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Instruction Manual

DDBioZ

DDBioZ

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Valid for the variants:

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Instruction Manual

1. Indications

DD Bio Z are dental mill blanks made from zirconium dioxide for the manufacture of dental prostheses. DD Bio Z can be used for all dental constructions with up to two bridge elements in the front and side tooth area.

2. Contraindications

Insufficient space, inadequate preparation, known intolerance to the contained components, insufficient oral hygiene.

3. Safety information

Please pay attention to the information in the safety data sheet in its current version. Avoid the inhalation of dust particles during processing. Wear protection gloves, safety glasses and a dust mask to prevent irritations of eyes, skin and respiratory system.

4. Handling and storage

The blanks are manufactured and tested to highest quality standards. Verify the integrity of the packaging and the blank before the first processing. Check if the content of the packaging corresponds to the declaration given on the label. Do not use damaged material. Storage only in the original containers in a cool and dry environment. Avoid vibrations, contaminations and contact with fluids.

5. Processing / Designing

Handling this medical device should only be allowed to trained staff. Monolithic restorations are suitable for bruxists – Generally, the following construction parameters need to be considered:

Indication	Schema	Wall thickness for all DD Bio Z variants (excl. ZS) [mm]	Wall thickness for DD Bio ZS Sirona inLab® [mm]	Connector cross section [mm ²]
Single crown	X	0.4	0.5	-
Telescope	X	0.5	0.5	-
Bridge 3-elements Anterior Posterior	XOX	0.5	0.5	> 7
		0.5	0.5	> 9
Bridge 4-elements Anterior Posterior	XOOX	0.6	0.7	> 9
		0.6	0.7	> 12

Legend: X = abutment tooth O = pontic

Please note: Depending on the construction, the connector cross section might need to be of larger dimensions. An oval connector cross section is desirable; the height of the connector is the decisive factor. Frameworks for ceramic veneering need to be designed in a way to support the veneering ceramic in the cusp region and allow a constant layer thickness. Shoulder or chamfer preparation is recommended.

6. Milling, sintering and further processing

Machining of blanks should always be conducted in the corresponding milling system. All information given by the manufacturer of the machine shall be taken into account. Information for Amann Girrbach users: Please find code (scale factor) for utilization of 98 mm discs in Amann Girrbach milling machines on blank margin. The frameworks need to be examined for visible defects. Do not process damaged or contaminated frameworks. The white framework can be colored before final sintering with DD Basic Shade / DD Pro Shade Z and DD Art Elements (pay attention to separate instruction manual). Pre-colored blanks are not allowed to use with patented coloring technologies.

Please pay attention to our separate sintering instruction!

Sintering cycle for normal furnace filling without cover:

- ↑ Heating up to 900°C (1652°F) with 8°C/min (46°F/min),
- dwell at 900°C (1652°F) for 30 minutes,
- ↑ heating up to final temp. 1450°C (2642°F) with 3°C/min (37°F/min),
- dwell at final sintering temperature 1450°C (2642°F) for 120 min,
- ↓ cooling to 200°C (392°F) with 10°C/min (50°F/min).

Avoid additional mechanical actions after sintering like blasting or grinding. If corrections are inevitable use water-cooled tools for the conditioning and make sure that there is no development of heat, which could lead to cracks in the material. Please work only with sharp, diamond-coated tools and very slight contact pressure. Areas that are under tension in clinical use (e.g. connectors) must not be processed. Do not separate in interdental regions. Avoid sharp edges.

7. Mounting

For anchoring we recommend conventional cementation with zinc oxide phosphates cements or glasionomer cements. Also luting composites may be used. Sufficient retentions and a stump with height of at least 3 mm are essential.

We reserve the right to make changes as a result of the continuous development of our products. Please find the current version of the instruction manual on our website: www.dentaldirekt.de This version replaces all previous versions.